

Product Data Sheet

Product Description: This cold kit formulation is optimized for the preparation of the small peptide radiopharmaceutical, ^{99m}Tc ---duramycin. The active ingredient in the kit is HYNIC-duramycin (Figure 1) in which the HYNIC chelator is attached to amino groups of duramycin. The conjugate contains 1 HYNIC molecule per duramycin molecule. Duramycin binds phosphatidylethanolamine (PE) at a 1:1 ratio with high affinity (K_d of 4-6 nM) and exclusive specificity.

Applications: ^{99m}Tc ---duramycin has been shown to be useful for *in vivo* SPECT imaging in animal models of sites where PE is exposed such as during apoptosis [Zhao, 2008]. In the cardiovascular area, ^{99m}Tc ---duramycin has been shown to be taken up in experimental models of myocardial infarction (MI) [Wang, 2015; Liu, 2016; Kawai, 2017], areas of atherosclerotic plaque [Liu, 2016; Hu, 2017; Kawai, 2017] and for imaging cardiotoxicity injury [Nakahara, 2017]. In the oncology area, it has been shown to allow early prediction of tumor response to therapy [Elvas, 2015 & 2016; Luo, 2016] therefore may be used as a screening tool for *in vivo* evaluation of new cancer drugs. It has also been used for whole body imaging to detect tissue injury from ionizing irradiation [Johnson, 2013] and may be useful to detect organ toxicity from known and experimental drug treatments [Kawai, 2017].

Product Formulation Procedure: For radiolabeling, ~30 mCi of ^{99m}Tc --pertechnetate in 0.5mL saline is added to the vial and excess vial pressure vented. The vial is then heated at 80°C in a lead-lined heating block for 30 mins and the radiolabel incorporation and radiochemical purity is assessed either by ITLC or radio-HPLC analysis [Elvas, 2015; Luo, 2016; Zhao, 2012]. Typical radiopurity of ^{99m}Tc -duramycin is >95%. Typical doses per rat and mouse are ~3--5 mCi and 1 mCi, respectively.

Blood Clearance Kinetics: The half-life for the fast clearance α -phase is ~3.6-4.1 mins whereas that of slow clearance β -phase is ~171-180 mins [Wang, 2015; Luo, 2016].

Size: 15 μg (~4-5 rat doses)

References:

- Zhao M et al., J. Nucl. Med. 49 (8), 1345-1352 (2008)
Wang L. et al., Nucl. Med. Biol. 42, 198-204 (2015)
Liu Z. et al. Nucl. Med. Biol., 43, 496-505 (2016)
Kawai H et al. in press (2017)
Hu Y et al. in press (2017)
Nakahara T et al. in press (2017)
Elvas F. et al., Mol. Imag. Biol. 17:838-847 (2015)
Elvas F et al., J. Nucl. Med. 57:805-811 (2016)
Luo R et al., Mol. Imaging, 15, 1-10 (2016)
Johnson SE et al. J. Nucl. Med. 54: 1397-1403 (2013)
Zhao M et al., Nucl. Med. Biol. 39:1006-1011 (2012)

Fig 1: Chemical structure of HYNIC-duramycin

