

## SUPERPARAMAGNETIC MAGNETITE

FeRex™ is a superparamagnetic magnetite. This reagent is classified as a darkening reagent acting through the T2 relaxation process. FeRex can be used as an MRI contrast reagent.

### References

1. *J Immunol Methods* 1982 **52**(3):353-367.
2. *FEBS Letters* 1984 **170**(2):282-238.
3. *J Aerosol Med Pulmonary Drug Delivery* 2008 **21**(4):335-341.

**NOT FOR HUMAN USE.**

### Superparamagnetic Magnetite

#### Catalog Number

MR-7200

FeRex .....\$ 175.00

10.0 ml sealed serum bottle at a concentration of 10 mg Fe/ml.  
NOT FOR HUMAN USE.  
Assay Method: MRI and EM.

