

SAMPLE VIALS FOR NEUTRON ACTIVATION PRODUCT LINE

BioPAL's polypropylene sample vials have been cleaned to remove sodium contamination as well as other trace elements and are calibrated for our automated activation and counting system. These vials are designed to contain biological samples, such as tissue and blood, during neutron activation. Customers **MUST** use these sample vials for BioPAL's neutron activation assay service.

Contaminant Free, Sample Vials

Customers **MUST** use BioPAL's sample vials for our neutron activation assay service. Customers are encouraged to contact BioPAL for more information concerning questions related to experimental protocol using neutron activation analysis.

Catalog
Number

V-0MS01

100 Polypropylene tissue-sample vials (2 ml) \$ 60.00 per pk

Because only 100 μ l sample of blood or urine is required, the 2 ml sample vials are ideal for GFR studies.

