

Applications and Functionalities

Catalog No.	Name	Application	Functionality
CL-30Q02-2	Molday ION	Blood pool agent	Hydroxyl
CL-30Q02-6	Molday ION (-)	Cell labeling with PLL	Hydroxyl
CL-50Q02-6A	Molday ION C6Amine	1. Cell labeling 2. Conjugation	Amino
CL-30Q02-6C	Molday ION Carboxyl Terminated	1. Cell labeling with PLL 2. Conjugation	Carboxyl
CL-50Q02-6A-50	Molday ION Rhodamine B	Cell labeling	Rhodamine B

Physical properties

Catalog No.	Name	~Size nm	~Zeta potential mV
CL-30Q02-2	Molday ION	30	-4.8
CL-30Q02-6	Molday ION (-)	30	-15
CL-50Q02-6A	Molday ION C6Amine	35	48
CL-30Q02-6C	Molday ION Carboxyl Terminated	35	-38
CL-50Q02-6A-50	Molday ION Rhodamine B	35	31

Volume, Concentration and Price

Catalog No.	Name	Volume ml	Concentration mg Fe/ml
CL-30Q02-2	Molday ION	2	10
CL-30Q02-6	Molday ION (-)	2	10
CL-50Q02-6A	Molday ION C6Amine	2	5
CL-30Q02-6C	Molday ION Carboxyl Terminated	2	2.5
CL-50Q02-6A-50	Molday ION Rhodamine B	2	2