



80 Webster Street  
Worcester, MA 01603 (USA)  
TEL: 508-770-1190

Web Site: [BioPAL.com](http://BioPAL.com)  
e-mail: [cs@BioPAL.com](mailto:cs@BioPAL.com)

## Product Data Sheet

---

Product Name:	Molday ION Carboxyl
Catalog Number:	CL-30Q02-7
Description:	Stable aqueous suspension of iron-based superparamagnetic nanoparticles (USPIO)
Application:	MRI contrast agents used primarily for vascular imaging
Packaging:	2 ml sealed serum bottle
Volume:	2.0 ml
Weight of Iron:	10 mg Fe/ml
Mean Particle Size:	20 nm
Core:	Magnetite
Mean Core Size:	8 nm
Coating:	Dextran
Number of Particles:	~2.2E17 per g
Density:	~1.25 g/cm <sup>3</sup>
Zeta Potential:	~ -20 mV
R1 (mM <sup>-1</sup> sec <sup>-1</sup> )	32.4
R2 (mM <sup>-1</sup> sec <sup>-1</sup> )	62.7
Storage Buffer:	dH <sub>2</sub> O
Autoclaved:	No
Long-Term Storage:	4 °C - <b>Do not freeze!</b>
Similar Products:	Ferumoxital, Feraheme
Expiration Date:	None

---